

High-Speed Overhauls



Standard overhaul includes replacement of bearings, O-rings and seals; concentricity check; cleaning of air and water lines; cleaning and polishing exterior; and replacement of end gasket. Testing under load at 35 to 40 psi.

	Warranty	Standard Overhaul
Midwest	6 mos	\$99
Star	6 mos	\$99
KaVo	6 mos	\$99
Lares	6 mos	\$99
NSK/Brasseler	6 mos	\$99
W&H/A-dec	6 mos	\$99
Bien-Air		Call for prices
All others		Call for prices

Autochuck handpieces with non-rebuildable chucks may require turbine replacement.

Replacement Turbines



No charge for turbine installation.

OE Turbines

	Warranty	Price
Midwest	6 mos	\$269
Star	1 year	\$325
Lares	6 mos	\$184
NSK/Brasseler	6 mos	\$199
Bien-Air		Call for prices
All others		Call for prices

Aftermarket Turbines

	Warranty	Price
Midwest	6 mos	\$169
Star	6 mos	\$169
KaVo	1 year	Starting at \$199
Lares		NA
NSK/Brasseler	6 mos	\$169
W&H/A-dec	6 mos	\$169
Bien-Air		NA
All others		Call for prices

Low-Speed Overhauls




Standard overhaul includes replacement of bearings, O-rings, seals and (depending upon make and model) all predetermined replacement parts required to achieve optimum performance. Cleaning and testing under load at 35 to 40 psi.

	Warranty	Standard Overhaul
Midwest	1 year	\$179
Star	1 year	\$179
KaVo	1 year	\$179
Lares	6 mos	\$179
NSK/Brasseler	6 mos	\$179
W&H/A-dec	6 mos	\$179
Bien-Air		Call for prices
All others		Call for prices

If the repair requires more than a standard overhaul, such as motor housing replacement or rotor housing, additional charges may apply.

Electric Overhauls



Product is overhauled and tested to meet new product specifications. The standard overhaul includes disassembly, a thorough cleaning and replacement of all necessary parts.

	Warranty	Standard Overhaul
KaVo		Free est.
NSK		Free est.
Star		Free est.
Bien-Air		Free est.
All others		Free estimates

Attachments



Overhaul includes replacement of bearings, gears, springs and seals (when applicable); cleaning; lubrication; and testing.

	Warranty	Standard Overhaul
Midwest Straight	6 mos	\$99
MW Contra Angle	6 mos	\$65
Star Straight	6 mos	\$99
Star Mtr-to-Angle	6 mos	\$75
KaVo		Starting at \$135
All others		Call for prices

Low-Speed Head Repair/Replacement



	Price
KaVo	\$135
Midwest Push Button (head overhaul)	\$75
Star Ball Bearing (head overhaul)	\$75
Generic Head Replacement	Starting at \$59

Non-rebuildable cartridges may require replacement. 3-6-month warranty.

Sonic Scalers



Sonic overhaul includes replacement of rotor, O-rings, seals and grommet; cleaning; clearing of air and water lines; and testing.

Star	\$159
All others	Call for prices

6-month warranty

Ultrasonic Scalers

Quick service (genuine Dentsply parts). \$50 declined estimate fee.

Most models including Dentsply & Parkell	Call for prices
--	-----------------

90-day warranty on parts and labor

Fiber Optic Replacement

	Price
Midwest Tradition, Quiet Air	\$159
Star 430	\$159
KaVo 640, 642, 635, 647, 649, 2307	\$240
NSK/Kinetic	\$240
All others	Call for prices

FACTORY-TRAINED

Bien-Air StarDental

LaRes
Your cutting edge

MIDWEST KAVO

HIGH-QUALITY REPAIRS

PROMPT TURNAROUND

FREE PICK-UP

LOCAL SERVICE

AUTOMATIC WARRANTY TRACKING