



Your Local Sharpening Service

Product Catalog & Reference Guide

Contact your LOCAL Sharpeez representative



Dear Customer,

Let me take a moment to introduce you to Sharpeez, a unique company that provides LOCAL hand instrument service. Sharpeez is dedicated to serving all your sharpening, retipping and new instrument needs.

Sharpeez is the sister company of Hayes Handpiece Company, a local handpiece repair company with 15 years experience. Like Hayes, Sharpeez has local offices to serve you—for better service and faster turnaround.

Here's what Sharpeez can do for you:

Sharpening



- We sharpen anything with a cutting edge including scalers, curettes, explorers, scissors and needle holders
- Sharper instruments make every procedure faster and easier, allowing you to maximize your patient load
- All major manufacturers

Retipping



- All 440 A U.S. stainless steel tips—the highest grade of steel for hand instruments
- All tips hardened to 55-57 on the Rockwell Hardness Scale so tips stay sharper for an extended period of time
- Keep your handle—if your instrument can be retipped, you can keep your existing handle
- All retipping guaranteed

New Instrument Sales



- High-quality instruments developed by Dental USA Inc.
- Handcrafted with the best material available
- Calibrated to the most precise specifications in the United States
- Choose from a variety of handles—high-quality steel to resist tarnishing, cracking and deformity

Please contact your LOCAL Sharpeez office for all your SHARPENING, RETIPPING or NEW INSTRUMENT needs!

Thanks for choosing Sharpeez!

Never a Dull Moment!



Suggestion Guide from Sharpeez

Best universal instrument for scaling & root planning:

3/4 Langer

Langer 3/4 = universal curette designed for maxillary arch

Langer 1/2 = universal curette designed for the mandibular arch

But both of the Langers can be used throughout the mouth where ever they adapt.

Anterior (centrals, laterals & canines):

Sub-0 , Gracey 1/2, 5/6, BH 1/2, YG 7/8, H6/H7, N5, N137

Posterior (molars & bicuspids):

YG 7/8, McCall 17/18, Nevi 2, 204S

Gracey 5/6, 7/8, 11/12 & 13/14, 17/18

Furcations:

Mini Graceys (60% of the time, furcation openings are too narrow to adapt a new curette, so mini's are a good choice.)

Barnhart 1/2

Colombia 4R/4L


Handle Options

	#	Descriptions	size
	1	Stainless Tapering Round	1/4"
	3	Stainless Solid Round	3/8"
	4	Stainless Ultra Lite Round	5/16"
	5	Stainless Octagon	1/4"
	6	Stainless Ultra Lite Round	3/8"
	7	Stainless Tapering Round	5/16"
	8	Stainless Short Lite Round	3/8"
	9	Stainless Short Octagon	3/8"
	10	Stainless Small Elevator Grip	1/2"
	11	Stainless Large Elevator Grip	1"
	12	Stainless Super ultra lite Satin Finish Round	3/8"

- We use the highest quality stainless steel for our handles
- Our handles are guaranteed against any tarnishing, cracking and deformity over time





Please inquire for further options

DIAGNOSTIC EXPLORERS



	HF	SHARPEEZ
S.E.	\$8.00	\$7.00
D.E.	\$11.25	\$10.00

DIAGNOSTIC PROBE-COLOR



	HF	SHARPEEZ
S.E.	\$15-19	\$14.00
D.E.	\$25.00	\$21.00

PERIODONTAL SICKLE/JACQUETTE SCALERS



Sharpeez
Your Local Sharpening Service
www.sharpeez.com

	HF	SHARPEEZ
D.E.	\$25.00	\$19.00

PERIODONTAL SICKLE SCALERS / CURETTES




Sharpeez
Your Local Sharpening Service
www.sharpeez.com

	HF	SHARPEEZ
S.E.		\$15.00
D.E.	\$25.00	\$19.00

PERIODONTAL SCALER & CURETTES

GOLDMAN FOX



CATTONI




COLUMBIA



ENGLISH PATTERN




UNIVERSAL BARNHART



	HF	SHARPEEZ
D.E.	\$25.00	\$19.00

PERIODONTAL UNIVERSAL CURETTES



LANGER (UNIVERSAL)



PROPHY GRACEY

SHANK LITTLE SHORTER THAN NOMAL GRACEY




Your Local Sharpening Service
www.sharpeez.com

	HF	SHARPEEZ
D.E.	\$25.00	\$19.00

PERIODONTAL


STANDARD GRACEY



AF-GRACEY (AVAILABLE SHORTER BLADE)



SUB-0 GRACEY



RIGID & SHORTER
BLADE
AVAILABLE

PROPHY GRACEY

SHANK LITTLE SHORTER THAN NOMAL GRACEY




Your Local Sharpening Service
www.sharpeez.com


	HF	SHARPEEZ
D.E.	\$23.00	\$19.00

PERIODONTAL


O'HEHIR CURETTES DEBRIDEMENT




HOES




FURCATION CURETTES




DARBY PERRY



FILES



FILES





Sharpeez
Your Local Sharpening Service
www.sharpeez.com


	SHARPEEZ FILES	SHARPEEZ CURETTES
D.E.	\$25.00	\$19.00

PERIODONTAL

SURGERY KNIVES








FILES



SHARPEEZ
\$24.00
\$29.00

PERIODONTAL SURGERY





Sharpeez
Your Local Sharpening Service
www.sharpeez.com


SHARPEEZ
\$21.00
\$22.00

ENDODONTIC PLUGGERS & CONDENSERS


RCP SET



ANTERIOR



POSTERIOR




	HF	SHARPEEZ
S.E.	\$17.00	\$15.00
D.E.	\$13.00	\$16.00

SURGICAL ELEVATORS



HEIDBRINK




SELDIN



Sharpeez
Your Local Sharpening Service
www.sharpeez.com

	HF	SHARPEEZ
S.E.	\$46.00	\$36.00


SURGICAL ELEVATORS



Sharpeez
 Your Local Sharpening Service
www.sharpeez.com

	HF	SHARPEEZ
S.E.	\$46.00	\$36.00


FORCEPS ENGLISH PATTERN






Your Local Sharpening Service
www.sharpeez.com

HF	SHARPEEZ
\$128.00	\$80.00

FORCEPS SPONGE & DRESSING




HAEMOSTATIC FORCEPS



Your Local Sharpening Service
www.sharpeez.com


HF	SHARPEEZ
\$47.00	\$27.00

CUTTERS & PLIERS









SHARPEEZ	
5601-5609	\$70.00
5610-5621	\$35.00

PLIERS



ORTHODONTIC INSTRUMENTS



Sharpeez
Your Local Sharpening Service
www.sharpeez.com

SHARPEEZ	
5601-5609	\$70.00
5610-5621	\$35.00

TISSUE, COTTON & DRESSING LIERS





Your Local Sharpening Service
www.sharpeez.com

ORDER FORM

Just fill out your order information and fax it to us or give to your local Sharpeez representative.

YOUR Office: _____

Contact: _____

Phone: _____

PAGE#	DESCRIPTION	PART NUMBER	SINGLE/DUAL END	QUANTITY
1	Diagnostic Explorers			
2	Diagnostic Probe Color			
3	Periodontal Sickle Jacquette Scalers			
4	Periodontal Sickle Scalers & Curettes			
5	Periodontal Scalers & Curettes			
6	Periodontal Universal Curettes			
7	Periodontal Gracey			
8	Periodontal Curettes/Hoes/Files			
9	Periodontal Knives			
10	Periodontal Surgical Curettes			
11	Endodonic Pluggers & Condensers			
12	Surgical Elevators			
13	Surgical Elevators			
14	Forceps English Pattern			
15	Forceps, Sponge & Dressing			
16	Haemostatic Forceps			
17	Cutters & Pliers			
18	Pliers			
19	Orthodontic Instruments			
20	Tissue, Cotton & Dressing Liers			

After filling out your order information fax it to us or give to your local Sharpeez representative.

